



new wide range -180°C to +110°C

EU 81 POSITION CRYOBOX

For Diagnostic & Clinical procedures

excellent performance and handling

- Universal fit for 0.5 ml up to 5 ml tubes
- “soft” PP for excellent handling
- Transparent Lid
- Stackable
- Available in 5 different colors
- Can handle tube heights from 1.6 to 6.7cm
- Can handle tube diameters of 12 mm
- Extra Robust & durable
- Can be use down to -180°C
- Liquid nitrogen resistance
- Can be use up to + 110°C
- Fits most freezer (cryo) storage systems
- Outer dimensions 13x13x4.5 cm



BIoplastics Cryo storage boxes allow visual examination of the box contents through the lid. The boxes are manufactured of durable polypropylene, which does not have the problem of becoming water saturated like cardboard boxes. The autoclavable, unbreakable design provides convenient storage for micro centrifuge tubes and cryo vials. The boxes can be stacked and have telescope lids. They are available in different colors. Working range from ambient temperatures up to 110 °C and down to -180°C. Accommodates all commonly used centrifuge and cryo tubes and fits most (Cryo) Storage Rack Systems.

Outer dimension 13.0 cm (L) x 13.0 cm (W) x 4.5 cm (H)


Cryo Storage Box 81 position

Outer dimension 13.0 cm (L) x 13.0 cm (W) x 4.5 cm (H)

Working range from ambient temperatures up to 110 °C and down to -180°C.

Orderno	Description	Package Size	Price
B10016	Cryo Storage Box 81 position	Natural, Box of 5	
B10017	Cryo Storage Box 81 position	Red, Box of 5	
B10018	Cryo Storage Box 81 position	Yellow, Box of 5	
B10020	Cryo Storage Box 81 position	Green, Box of 5	
B10021	Cryo Storage Box 81 position	Blue, Box of 5	





MADE IN THE NETHERLANDS

Products are designed and manufactured in the Netherlands

Usage range -180°C to +110°C

Head Quarters:

BIOplastics BV
Rötscherweg 61
Entrance D
6374 XW Landgraaf
The Netherlands

Phone: +31 (0)45 533 87 50
Fax: +31 (0)45 533 87 96
Website: www.bioplastics.com
E-mail: info@bioplastics.com

Subsidiary:

BPCTi Inc.
2933 South Miami Blvd.
Suite 121
Durham, NC 27703, USA

Phone: 919 806 8811
Fax: 919 806 2014
E-mail: info@bpcti.com
Website: www.bpcti.com

Your local distributor: