


single 0.1 ml (q)PCR tubes and 8-tube-strip with single attached caps which fits Roche 480*, ABI Fast Cyclers, 7500, Veriti, StepOne StepOnePlus, Bio-Rad CFX and many other cyclers.

For Diagnostic & Clinical procedures


- Fits ABI Fast Cyclers, Roche 480, Bio-Rad CFX and many other qPCR cyclers(see table next page)
- Designed to perform diagnostic qPCR and PCR reactions
- Available in single tubes and 8-tube strip with single attached caps
- Attached caps prevent cross contamination
- Wide area indented Optical Cap enables maximum signals and minimizes optical area contamination
- Light frosted outer tube-surface enhances qPCR signals
- Extra marking areas and auto orientation holes (8-tube strip) minimizes identification mistakes
- Extreme Uniform wells lead to reproducible results
- Low evaporation due to clever design (can be used with 5 ul)



light frosted surface enables signal enhancement in qPCR



Indented Optical Cap prevent area contamination



orientation & "click" hole

EU 0.1 ml Thin-wall tube with optical indented wide area cap. Low Profile (0.1 ml)

Fits almost all cycler models and qPCR (fast) cycler models which accept low profile products. (ABI Veriti(fast), Roche 480*, Bio-Rad, Eppendorf)(*requires specific adaptor) Can be used in most 0.2 ml Fast and Real-Time thermal cyclers .

Orderno	Description	Package Size	Price
B77201	EU 0.2 ml Thin-wall tube, Low Profile, natural,	bag, 1000	check our website or contact your local distributor

8-tube strip with single attached optical wide area indented cap. Low Profile (0.1 ml)

Orderno	Description	Package Size	Price
B72811**	EU 0.2 ml Thin-wall 8-tube strip with single attached optical caps, LP, naturalbag, 120		check our website or contact your local distributor

**Not available in the USA, alternatives available: B77001/B79501





(q)PCR Cycler reference chart for 0.1 ml tubes with optical indented wide area cap, Low Profile (0.1 ml)
 or use the interactive search engine cycler to product fit at: www.bioplastics.com/services/filter.aspx


All Cycler brands with height adjustable lids: **YES**



	Mini Opticon	GeneCycler	iCycler	MyCycler	PTC 100	PTC 200	PTC 220	PTC 221	PTC 225	Chromo-4	iCycler iQ	IQ 5	MyiQ	CFX 96	C1000	S1000
Applied Biosystems (fast)						✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Bio-Rad	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	MasterCycli.	Master Grad.	Master ep	ep Realplex												
Eppendorf	✓	✓	✓	✓												
	480 96 block with adaptor B79480															
Roche (requires adaptor B79480)	✓															



Roche 480 adaptor B79480



Products are designed and manufactured in the Netherlands. SAMPLES on request and/or purchase with money back guarantee* (* if proven to be not working) Contact your local representative.

All rights reserved. 052010V1

Trademarks; owner marks; > Applied Biosystems: - Veriti - StepOne - StepOnePlus ; Bio-Rad: - Opticon - iCycler - MyCycler - Chromo-4 - C1000 - S1000 ; Eppendorf: - MasterCycler - RealPlex ep ; Roche: -Roche 480

Head Quarters:

(R&D and Manufacturing)
 BIOplastics BV
 Rötischerweg 61
 Entrance D
 6374 XW Landgraaf
 The Netherlands

Phone: +31 (0)45 533 87 50
 Fax: +31 (0)45 533 87 96
 Website: www.bioplastics.com
 E-mail: info@bioplastics.com

Subsidiary:

BPCTi Inc.
 2933 South Miami Blvd.
 Suite 121
 Durham, NC 27703, USA

Phone: 919 806 8811
 Fax: 919 806 2014
 E-mail: info@bpcti.com
 Website: www.bpcti.com

Your local distributor:

Please contact
 BIOplastics Headquarters
 The Netherlands
 if distributor informaton is lacking