


BIOplastics

48-well (q)PCR plates

Easy to cut



B71601



B57651



B57651

Availability

- 0.1 ml low profile, semi skirted
- 0.2 ml regular profile, non skirted
- Frosted & white, laser marked
- Light frosted, white & ultra clear
 For resp. high/highest S/N ratios in qPCR use light frosted/white

In short

- Semi skirted version can be cut per three rows offering flexible usage
- · Sample identification
- · Laser marking possibilities
- All 8-cap strips fit
- White plates maximize S/N ratios & avoid signal contamination from the block

Compatible cyclers

0.1 ml: Fits ABI StepOne, Veriti-Fast, Bio-Rad CFX, Chromo-4, (Mini) Opticon, Eppendorf Real Plex and many other PCR and qPCR cyclers. 0.2 ml: Fits almost all (q)PCR cyclers that accept regular profile products.

See compatibility list at www.bioplastics.com.

(q)PCR Perfect

Highly reproducible results
Low evaporation
Low binding
Stackable
No nesting
Robust
Pyrogen, DNa(se) & RNa(se) free

Products are designed and manufactured in The Netherlands. Samples on request.

© 2019 BIOplastics B.V. All rights reserved

www.bioplastics.com +31(0)455338750 order@bioplastics.com




BIOplastics


Ordering information

48-well (q)PCR plates

| | Light frosted | White | Ultra Clear | Packaging |
|--------------------------------|---------------|---------------|---------------|------------|
| 48 x 0.1 ml semi skirted plate | <u>B71601</u> | <u>B71609</u> | | 50 plates |
| 48 x 0.2 ml non skirted plate | | <u>B71509</u> | <u>B71501</u> | 50 plates |
| | | | | |
| Closure options | | | | |
| EU optical wide 8-cap strip | <u>B57801</u> | | | 120 strips |
| Tear Off 8-cap strip mat | <u>B57651</u> | | | 50 mats |







B57801



B71609

Products are designed and manufactured in The Netherlands. Samples on request.

© 2019 BIOplastics B.V. All rights reserved