

Budget friendly

# Most universal (q)PCR 96 well plate


Fits most PCR and qPCR Cyclers

For Diagnostic & Clinical procedures

- Designed to perform Diagnostic qPCR and PCR reactions in volumes down to 5 µl
- Fits ABI, Bio-Rad, MJ Research, Eppendorf and many other qPCR and regular cyclers (see table on page 2)



FEATURES	ADVANTAGES	BENEFITS
Frosted plates	Higher signal/noise ratio's	Lower copy detection
Special blend of polypropylene	Anti-static tubes	Maximum sample recovery
	Non binding for biomolecules	Low copy detection
	Flexible non cracking tubes	No loss of sample
Uniform wall thickness	Uniform heating	Consistent reaction conditions
	No risk of cracking	No loss of sample
	Minimize evaporation	Consistent reaction conditions
Frosted writing areas	Easy sample marking	Traceability of samples
Optimal tube-block contact angle	Better heat transfer	Optimal reaction conditions
Thin wall lower part tube	Maximal heat transfer	Optimal reaction conditions
Thick wall upper part tube	Minimal evaporation	Consistent reaction conditions
Numeric column markings	Easy sample location	Traceability of samples
Alphabetic row markings	Easy sample location	Traceability of samples
Vertical orientation lines (5)	Easy sample location	Traceability of samples
Horizontal orientation lines (3)	Easy sample location	Traceability of samples



light frosted surface enables signal enhancement in qPCR

Orientation Lines

## € 42.50 for 25 plates

### EU Semi-skirted light frosted 96 x 0.2 ml (q)PCR thin-wall plate. Regular Profile

These EU 96 well regular profile plates are semi-skirted and designed for (q)PCR applications. Closure can be accomplished with either EU (optical), strip caps or seals.

Order no	Description	Package Size	Price €
AB19700	EU 96 x 0.2 ml Thin-wall plate Semi-skirted, RP, LF, natural	.25 plates	42.50

### Closure options

Order no	Description	Package Size	Price €
B57801	EU Optical Wide 8 cap-strip, with indented flat cap, natural	.bag, 120	18.70
B57801B	EU Optical Wide 8 cap-strip, with indented flat cap, natural	.bag, 300	47.60
157300	Opti-Seal, optical, disposable, adhesive	.100 sheets	91.62

For other closure models and cap-strip options see [www.bioplastics.com](http://www.bioplastics.com)



(q) PCR Cycler Compatibility List for AB19700 plates		Other Cyclers	
Most popular cyclers		Other Cyclers	
Applied Biosystems	310	Abbott	m2000rt
Applied Biosystems	3100	AcuGen Systems	AG-9600
Applied Biosystems	3130	Analytik-Jena	FlexCycler
Applied Biosystems	3700	Apollo	ATC201
Applied Biosystems	3730	Apollo	ATC401
Applied Biosystems	7300 Real-Time PCR System	Apollo	ATC701
Applied Biosystems	7500 Real-Time PCR System	Barnstead Thermolyne	Amplitrion I
Applied Biosystems	7900HT Real-Time PCR System	Barnstead Thermolyne	Amplitrion II
Applied Biosystems	ABI PRISM 7000 SDS	Beckman	CEQ 8000
Applied Biosystems	ABI PRISM 7700 SDS	Beckman	CEQ 8800
Applied Biosystems	GeneAmp PCR System 2700	Bioer	Life Express TC-96/T/H(a)
Applied Biosystems	GeneAmp PCR System 2720	Bioer	XPcyler TC-XP-A
Applied Biosystems	GeneAmp PCR System 5700	Bioer	XPcyler TC-XP-G
Applied Biosystems	GeneAmp PCR System 9600	Bioneer	Exicycler 96
Applied Biosystems	GeneAmp PCR System 9700	Bioneer	MyGenie 96
Applied Biosystems	GeneAmp PCR System 9700 Aluminum	Boeco	TC-PRO 96G
Applied Biosystems	GeneAmp PCR System 9700 Gold	Corbett Research	FTS 96
Applied Biosystems	Solid3	Corbett Research	FTS A 96
Applied Biosystems	Veriti™ 96-Well Thermal Cycler	Corbett Research	Palm-Cycler
Bio-Rad Laboratories	C1000	Corbett Research	PC 96
Bio-Rad Laboratories	CFX96 Real-Time PCR Detection System	DuPont Qualicon	BAX System Q7
Bio-Rad Laboratories	DNA Engine 96v	eENZYMÉ	MyGene MG-96
Bio-Rad Laboratories	DNA Engine Chromo4	Ericom	Deltacycler II
Bio-Rad Laboratories	Dyad 96v	Esco	Swift™ Maxi
Bio-Rad Laboratories	Dyad Chromo4	Esco	Swift™ MaxPro
Bio-Rad Laboratories	Dyad Disciple 96v	Esco	Swift™ Spectrum 96
Bio-Rad Laboratories	Dyad Disciple Chromo4	EuroClone	One Advanced
Bio-Rad Laboratories	Tetrad 2 96v	EuroClone	One Gradient
Bio-Rad Laboratories	Tetrad 2 Chromo4	EuroClone	peqSTAR 96 HPL
Bio-Rad Laboratories	iCycler	EuroClone	peqSTAR 96 ML
Bio-Rad Laboratories	iQ Real-Time PCR	EuroClone	peqSTAR 96 Universal
Bio-Rad Laboratories	iQ5 Real-Time PCR	Finnzymes Instruments	TCA0096 Arktik
Bio-Rad Laboratories	MyCycler	Gene Technologies	G-Storm GS1/GS2/GS3/GS4
Bio-Rad Laboratories	MyiQ Real-Time PCR	Helena Biosciences	Phoenix Thermal Cycler
Bio-Rad Laboratories	PTC-100	INHECO	Heater Olly 96
Bio-Rad Laboratories	S1000	Lab-Line Instruments	Programmable Thermal Block
Eppendorf	Mastercycler	Labnet	MultiGene™
Eppendorf	Mastercycler ep gradient	Landgraf	Varius M
Eppendorf	Mastercycler ep gradient Silver	LongGene	Gradient PCR
Eppendorf	Mastercycler ep realplex	LongGene	L Series 96+
Eppendorf	Mastercycler ep realplex Silver	LongGene	Multi-Purpose PCR
Eppendorf	Mastercycler pro	LongGene	MyGene PCR
Eppendorf	Mastercycler pro S	MWG-Biotech	Primus 96 Advanced
Eppendorf	Thermomixer comfort	MWG-Biotech	Primus 96 Plus
MJ Research Inc.	PTC-100	MWG-Biotech	Primus Dualblock PCR-System
MJ Research Inc.	PTC-100 AgV	MWG-Biotech	Primus Multiblock PCR-System
MJ Research Inc.	PTC-200 DNA Engine 96v	MWG-Biotech	Primus PCR-System
MJ Research Inc.	PTC-200 DNA Engine Chromo4	Peqlab	peqSTAR 96 *
MJ Research Inc.	PTC-220 DNA Engine Dyad 96v	Peqlab	Primus 96 Advanced
MJ Research Inc.	PTC-220 DNA Engine Dyad Chromo4	ProGen	Thermal Cycler
MJ Research Inc.	PTC-221 DNA Engine Disciple 96v	PTI	FluoDia T70
MJ Research Inc.	PTC-225 DNA Engine Tetrad 96v	Scinics	EZ-Cycler MP-96
MJ Research Inc.	PTC-225 DNA Engine Tetrad Chromo4	SensoQuest	Labcyler Themoblock 96
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2 96v	Techne	Flexigene
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2Chromo4	Techne	Gene-E/Genius
Stratagene	Mx3000P™, Mx3005P™	Techne	TC-412/512/PLUS
Whatman Biometra	T1 Thermocycler + Combi	Techne	Touchgene Gradient
Whatman Biometra	TGradient	Thermo Hybaid	MultiBlock System
Whatman Biometra	TOptical	Thermo Hybaid	Omn-E
Whatman Biometra	TProfessional	Thermo Hybaid	Omnigene TR3 CM220
Whatman Biometra	TProfessional Basic + Basic Gradient	Thermo Hybaid	PCR Express, Px2, PxE
Whatman Biometra	TRobot	Thermo Hybaid	TouchDown
Whatman Biometra	UNO II	Transgenomics	Wave MD System 4000
Whatman Biometra	UNO-thermoblock	Xian Tianlong Science	DTC-4, D49

**Trademarks**

(registered) Trademarks and Owner: ABI Prism 5700, 700, 7000, 7700 and 7900 owner The Perkin Elmer Corp; Abbot, M2000rt owner Abbot Corp.; Opti-Seal BIOplastics (logo), EU, BIOplastics owner BIOplastics BV; Analytik-Jena, FlexCycler owner Analytik-Jena; CEQ owner Beckman; Corbett Research, PalmCycler owner Corbett Research; Dupont, BAQ Q7 owner Dupont Corp.; Eppendorf, MasterCycler owner Eppendorf AG; LightCycler owner Roche; MegaBACE owner GE Healthcare; Mx3000, Mx3000P, Mx3005P, Mx4000 owner Stratagene, Opticon, iCycler iQ, iCycler, iQ5 owner Bio-Rad laboratories Inc.; Quantica owner Techne Ltd; StepOne, Veriti owner Applera Corp. or it's subsidiaries; Thermo Hybaid owner Thermo Corp.; Wave owner Transgenomics

Products are designed and manufactured in the Netherlands. SAMPLES on request and/or purchase with money back guarantee\* (\* if proven to be not working) Contact your local representative. All rights reserved. 11070V1 ©BIOplastics BV

**Headquarters:**

(R&D and Manufacturing)  
 BIOplastics BV  
 Rötscherweg 61  
 Entrance D  
 6374 XW Landgraaf  
 The Netherlands

Phone: +31 (0)45 533 87 50  
 Fax: +31 (0)45 533 87 96  
 Website: www.bioplastics.com  
 E-mail: info@bioplastics.com

**Subsidiary:**

BPCTi Inc.  
 2933 South Miami Blvd.  
 Suite 121  
 Durham, NC 27703, USA

Phone: 919 806 8811  
 Fax: 919 806 2014  
 Website: www.bpcti.com  
 E-mail: info@bpcti.com

