


# EU Thin Wall (q)PCR 24 well semi-skirted plate


## For Diagnostic & Clinical procedures

e x c e l l e n t p e r f o r m a n c e a n d h a n d l i n g

- Perfect fit to nearly ALL PCR & qPCR cyclers  
(see table next page)
- Fits ABI 7300 and 7500 Real Time Cyclers
- Designed to perform qPCR and PCR reactions
- Cutable
- Low evaporation due to clever design
- Easy labeling and orientation, can be color coded by “clicking” the grid onto the plate(s)
- Extreme Uniform wells lead to reproducible results
- Can be closed using any of the BIOplastics cap strips
- Can be placed, transported & stored in Tube Support grid, Station Support Insert or Work Rack S-96



light frosted surface enables signal enhancement in qPCR



orientation & “click” hole

Cut positions

### EU Semi Skirted Thin-wall 24 x 0.2 ml plate, 3 x 8

Can be used in most 0.2 ml block regular and Real-Time thermal cyclers  
(see table next page )

**Orderno Description**

B50240	24 x 0.2ml Plate, Natural, Box of 100 plates, Semi-skirted Thin-wall
B57801	EU Wide Optical 8-Cap Strip, Wide Indented Flat Cap
B75701	EU Optical 8-Cap strip
B69351	EU Pre-Post Tube Support Grid Wide, Natural, Re-useable

**Package Size**

Box, 100 plates
Natural, Bag of 120 Strips
Natural, Bag of 120 Strips
Natural, 8 Grids


**Price**

(q)PCR Cycler reference chart for Breakable Semi Skirted Plate

All Cycler brands with height adjustable lids: **YES**



		310	3100	3130	3700	3730	2400	9600	2700	2720	7000	7300	7500	7700	7900	9700	Veriti™ 96
<b>Applied Biosystems (non fast)</b>		X	X	X	X	X	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
		GeneCycler	iCycler	MyCycler	PTC 100	PTC 200	PTC 220	PTC 221	PTC 225	Chromo-4	iCycler IQ	IQ 5	MyIQ	CFX 96	C1000	S1000	
<b>Biorad</b>		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
		MasterCycl.	Master Grad.	Master ep	ep Realplex												
<b>Eppendorf</b>		✓	✓	✓	✓												
		Mx 3000P	Mx 3005P														
<b>Stratagene</b>		✓	✓														



Products are designed and manufactured in the Netherlands. All rights reserved. 022009V1

Trademarks; owner: marks -> Applied Biosystems: - ABI Prism - Veriti - StepOne - StepOnePlus ; Bio-Rad: - Oplicon - iCycler - MyCycler - Chromo-4 - C1000 - S1000 ; Eppendorf: - MasterCycler - RealPlex ep ; Stratagene: - Mx3000P - Mx3005P

**Head Quarters:**

(design and Manufacturing)

BIoplastics BV

Rötscherweg 61

Entrance D

6374 XW Landgraaf

The Netherlands

Phone: +31 (0)45 533 87 50

Fax: +31 (0)45 533 87 96

Website: www.bioplastics.com

E-mail: info@bioplastics.com

**Subsidiary:**

BPCTi Inc.

2933 South Miami Blvd.

Suite 121

Durham, NC 27703, USA

Phone: 919 806 8811

Fax: 919 806 2014

E-mail: info@bpcti.com

Website: www.bpcti.com

**Your local distributor:**

