


4-tube strip

For Rotor-Gene® systems

Consistant Ct values



Availability

- 0.1 ml low profile
- Clear transparant

In short

- Individual attached caps minimize cross contamination
- Optimal tube-fit to rotor for optimal reaction conditions
- Volumes down to 15 µl
- Identification possibilities
- Fully traceable by 2D laser marked barcode



Compatible cyclers

All Rotor-Gene® (q)PCR cyclers from QIAGEN

(q)PCR Perfect

Highly reproducible results
 Low evaporation
 Low binding
 Stackable
 No nesting
 Robust
 Pyrogen, DNase & RNase free

Ordering information

4-tube strip

		Package size
EU 0.1 ml optical 4-tube strip	B84001	250 strips, 1000 tubes
EU 0.1 ml optical 4-tube strip 2D laser marked barcode	B85001L	250 strips, 1000 tubes

Designed and manufactured in The Netherlands. Samples on request.

© 2019 BIOplastics B.V. All rights reserved / 110719

www.bioplastics.com

+31(0)455338750

order@bioplastics.com