


Tear Off products

96-well Tear Off plate mat




Availability

- 0.1 ml low profile
 - 0.2 ml regular profile
 - 12 x 8-tube strip mat
 - 8 x 12-tube strip mat
 - Light frosted, white and clear
- For qPCR applications light frosted or white products are recommended

In short

- Create different well and strip formats
 - Easy to tear off
 - Tube strips are barely attached to each other
 - Strips don't bend or break
 - Easy identification of samples and strip orientation
 - Laser marking possibilities
 - Compatible with shell frame grids
- Tube mats can be clicked in specific (skirted) grids to create a plate assembly for specific (q)PCR instruments and/or robotic applications



Compatible cyclers

0.1 ml: fits all fast 0.1 ml (q)PCR thermocyclers and all 0.2 ml thermocyclers with adjustable lid
 0.2 ml: fits all fast 0.2 ml (q)PCR thermocyclers
 Compatible cycler list can be found on www.bioplastics.com



Closure

96-well Tear Off plate mats can be closed with any of the BIOplastics 8- or 12-cap strips, cap mats and Tear Off cap strip mats.

(q)PCR Perfect

Highly reproducible results
 Low target copy detection
 Low evaporation, low binding
 Stackable, no nesting
 Robust
 Pyrogen, DNase & RNase free

Designed and manufactured in The Netherlands. Samples on request.

© 2019 BIOplastics B.V. All rights reserved / 110219

www.bioplastics.com

+31(0)455338750




order@bioplastics.com

Ordering information

96-well Tear Off products

Product description	Low profile 0.1 ml	Regular profile 0.2 ml	Packaging
8 x 12 tube strip mat light frosted	B59151	B58151	25 plates/200 strips
8 x 12 tube strip mat white	B59159	B58159	25 plates/200 strips
8 x 12 tube strip mat ultra clear*	K59151	K58151	25 plates/200 strips
Optical Tear Off 12-cap strip mat (8 x 12)	B58651	B58651	50 mats/400 strips
12 x 8 tube strip mat light frosted	B59001-1	B58001-1	25 plates/300 strips
12 x 8 tube strip mat white	B59009-1	B58009-1	25 plates/300 strips
12 x 8 tube strip mat ultra clear*	K59001-1	K58001-1	25 plates/300 strips
Optical Tear Off 8-cap strip mat (12 x 8)	B57651	B57651	50 mats/600 strips

*Available in five colors



Designed and manufactured in The Netherlands. Samples on request.

© 2019 BIOplastics B.V. All rights reserved / 110219