

BIOplastics

(q)PCR Closure options

8-cap strips | 12-cap strips | 96-cap strip mats | seals

All cap strips

- Can close all BIOplastics (q)PCR tubes, strips & plates
- Are available in ultra clear

8-cap strips



B75701-1 Extra robust, indented flat cap-strip Great solution for storing samples



B79501
8-single attachable indented caps Ideal for real-time applications in combination with white or frosted tube strips
Can close BIOplastics products: B77001-11
K77001-11 B 77101-11
K77101-11



C79701 Semi domed cap-strip

8 cap-strips with

Optical wide area indented caps

Indented caps greatly reduce potential contamination by hands



B79701-1 Robust, yet flexible cap-strip



B79721
Flexible strip through thin cap to cap connections
Easy to cut with scissors or knife



Strip with flexible connections between the caps
Extra secure closure



B57801 Extra robust strip

Designed and manufactured in The Netherlands. Samples on request.

© 2019 Bioplastics B.V. All rights reserved / 280619

www.bioplastics.com

+31(0)455338750

order@bioplastics.com



BIOplastics

(q)PCR Closure options

8-cap strips | 12-cap strips | 96-cap strip mats | seals

12-cap strips



B57861
Optical wide area indented flat cap-strip Robust, yet flexible



B57821 Optical wide area indented flat cap-strip Robust strip



B56501 Robust, indented flat cap-strip Great solution for storing samples


Cap-plates & mats



96-cap Transformer ™ plate
Easy to cut into any required format



96-cap plate with wide area indented flat caps
Easy to cut with scissors into any required format



B57751
Optical cut-off 8-cap strip mat (12 x 8)
Strips are attached to each other with hinges
Can be cut in one or multiples 8-cap strips or two or multiples 12-cap strips



B57651 12 x 8-cap
B58671 8 x 12-cap
Tear Off strip mat
Strips are slightly attached to each other so
you can easily tear off one
or multiples 8- or 12-cap
strips

Designed and manufactured in The Netherlands. Samples on request.

© 2019 Bioplastics B.V. All rights reserved / 280619

www.bioplastics.com +31(0)455338750 order@bioplastics.com



BIOplastics

Ordering information

(q)PCR Closure options

Seals



Optical seals are very easy to apply & remove Three variations:

157300 Optical sealing sheet Optical sealing sheet

<u>157400</u>

Optical sealing sheet Laser Marked Optical sealing sheet | devidable in

strips to close (multiples of) 8 wells

Product description	Order number	Packaging
8-cap strips		
EU robust indented flat 8-cap strip	<u>B75701-1</u>	120 strips
EU 8-single attachable indented caps	B79501	120 strips
EU semi domed 8-cap strip*	<u>C79701</u>	120 strips
EU optical robust flat 8-cap strip	B79701-1	120 strips
EU optical flat 8-cap strip	B79721	120 strips
EU optical wide area 8-cap strip	<u>B57881</u>	120 strips
EU optical wide area 8-cap strip	<u>B57801</u>	120 strips
12-cap strips		
EU optical standard flat 12-cap strip	B57861	80 strips
EU optical wide area 12-cap strip	B57821	80 strips
EU optical indented flat 12-cap strip*	<u>B56501</u>	80 strips
mats/plates		
EU optical wide area 96-cap transformer plate	<u>B57701</u>	25 plates
EU optical side area 96-cap plate	B57601	25 plates
Optical cut off 8-cap strip mat (12 x 8)	<u>B57751</u>	50 mats/600 strips
Optical Tear Off 12-cap strip mat (8 x 12)	B58671	50 mats/400 strips
Optical Tear Off 8-cap strip mat (12 x 8)	<u>B57651</u>	50 mats/600 strips
seals		
Optical sealing sheet	<u>157300</u>	100 sheets
Optical sealing sheet, laser marked	<u>157300L</u>	100 sheets
Devidable Optical sealing sheet	<u>157400</u>	100 sheets

^{*}Available in different colors

Designed and manufactured in The Netherlands. Samples on request.

© 2019 Bioplastics B.V. All rights reserved / 280619

www.bioplastics.com +31(0)455338750

order@bioplastics.com